SUPLASYN[®] m.d.

STERILE SODIUM HYALURONATE SOLUTION Viscoelastic Supplement for Small Synovial Joints

HYALURONIC ACID (HA) is a normal component of synovial fluid and plays a central role in maintaining the physiological internal environment of the joint.

COMPOSITION: Viscoelastic solution of a defined molecular fraction of purified hyaluronic acid produced by fermentation. (*m.d.* stands for mini-dose).

Each syringe contains:

Hyaluronic acid sodium salt 7 mg
Excipients q.s 0.7 ml

PROPERTIES: Hyaluronic acid is the prototype of a wide range of saccharide biopolymers (glycosaminoglycans or mucopolysaccharides), important components of all extracellular tissue structures, including cartilage and synovial fluid. The active substance of SUPLASYN[®] *m.d.* is a specific hyaluronic acid fraction of a defined molecular chain length with a high degree of purity. The viscoelastic properties of SUPLASYN[®] *m.d.* are similar to synovial fluid.

INDICATIONS: As a supplement for small synovial joints. By improving rheological properties of the synovial fluid, SUPLASYN[®] *m.d.* has been shown to be beneficial in osteoarthritis for the management of pain and the improvement in the physical function of joints. For larger joints, SUPLASYN[®] syringes of 20 mg / 2 ml should be used.

CONTRAINDICATIONS: Do not administer to patients with a history of allergy to the components. No other contraindications to the intra-articular use of SUPLASYN[®] *m.d.* are known at this time.

DOSAGE AND ADMINISTRATION: Depending upon joint size, up to 0.7 ml may be administered intra-articularly. The recommended schedule is once a week for three weeks, but supplemental injections can be administered depending on the patient's clinical condition. Using aseptic technique, introduce up to 0.7 ml intra-articularly into the affected joint. More than one joint may be treated at the same time. The syringe contents are for single use only. Discard any unused portion of the syringe. To use the pre-filled syringe, remove the Luer lock cap, attach a suitable cannula (for example a 25 G) and secure it by turning slightly.

PRECAUTIONS: Transient short duration pain may occur following intra-articular introduction. The affected joint may show a mild reaction such as pain, feeling of heat, redness or swelling. Symptoms subside within days for most patients. In rare cases, stronger local reactions and/or fever may occur. Discontinue use if adverse reactions are experienced. As no clinical evidence is available on the use of hyaluronic acid in children, pregnant or lactating women, treatment with SUPLASYN[®] *m.d.* is not recommended in these patients. Avoid using SUPLASYN[®] *m.d.* with instruments sterilised with quaternary ammonium salts solutions.

WARNING: KEEP OUT OF REACH OF CHILDREN. DO NOT USE IF BLISTER IS DAMAGED. TO BE USED BY A PHYSICIAN ONLY.

PACKAGING: Available 0.7 ml syringes, packaged individually.

STORAGE: Store between 4°C and 30°C. DO NOT FREEZE. BRING TO ROOM TEMPERATURE BEFORE INJECTION.

Date of preparation: December 2003.

Manufactured by:
© Bioniche Pharma Group Ltd.
Inverin, County Galway
Republic of Ireland
www.suplasyn.com

Distributed by:
Bioniche Pharma (Canada) Limited
Belleville, Ontario K8N 5J2
1-888-258-4199

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